

Analytical Sensitivity and Specificity of IMMUNOCOMB® Antibody Test Kits

Product		Sensitivity	Specificity	Study site
 Canine Parvovirus & Distemper virus IgM	Parvovirus	91%	91%	Virology Laboratory, College of Veterinary Medicine. Auburn University, AL, USA
	Distemper virus	94%	96%	
 Canine Ehrlichia canis		100%	94%	Hebrew University of Jerusalem, Israel
 Canine Leptospira		98%	97%	Catholic University of Temuco, Temuco, Chile
 Canine Brucella		98%	93%	Biogal Galed Labs
 Feline Toxoplasma gondi & Chlamydia psittaci	T.gondi	92%	97.37%	Agrolabo Spa, Italy
	C.psittaci	94.74%	99%	
 Feline Corona Virus / FIP.		100%	100%	University of Glasgow ,UK